


E.O.S. DEV., FIX AND FIXER

Universal, environmental friendly developer and fixers



- Designed to be used for automatic processing in conventional systems
- Combination of E.O.S. Fix, E.O.S. Fixer and E.O.S. Dev. enables on-line silver recovery
- Screw caps with built-in cutting edges for easy opening
- Bottles fit for automatic mixing units



E.O.S. DEV. IS AGFA HEALTHCARE'S UNIVERSAL ECOLOGICAL HARDENING DEVELOPER. E.O.S. FIX / E.O.S. FIXER IS AGFA HEALTHCARE'S UNIVERSAL LOW-ODOR FIXER.

THE CHEMICALS ARE SPECIALLY DESIGNED TO BE USED WITH AGFA HEALTHCARE'S ECOLOGICALLY OPTIMISED SYSTEM (E.O.S.).

E.O.S. Dev. for automatic processing

E.O.S. Dev. is Agfa HealthCare's universal ecological hardening developer for automatic processing in conventional systems.

E.O.S. Dev. yields both outstanding sensitometric and physical properties.

Depending on the film type used and the throughput -average processed square metres per day- replenishment may vary from 300 to 600 ml/m².

E.O.S. Fix and E.O.S. Fixer, universal low-odor fixers

Agfa HealthCare's E.O.S. Fix and E.O.S. Fixer are ecologically optimized two-part hardening fixers. Especially designed to be used in combination with E.O.S. Dev.

Depending on the mixer type, E.O.S. Fix (2 x 25 l) or E.O.S. Fixer (2 x 20 l) is to be used.

Ecologically optimized

E.O.S. Dev. is an ecologically optimized derivate of the G138i. Together with the E.O.S. Fix or the E.O.S. Fixer, it offers you true ecological soundness. On-line silver recovery is no problem, and guarantees you to be able to comply with all current and pending environmental legislations.

Ecologically improved packaging

E.O.S. Dev., E.O.S. Fix and E.O.S. Fixer are shipped in cartons made of 60 % recycled cardboard. Two handles at both sides of the carton facilitate handling and transport. The concentrated liquids are packed in translucent polyethylene (PE) bottles with PE ECO seals. The screw caps and the stripe at the top of the label are colored for identification purposes. Developer bottles have a red cap, except for part C, which has a transparent cap to distinguish it from part B. The fixer bottles have blue caps. The colored components, however, have been reduced as much as possible and the bottles show no aluminium remnants after use. Both aspects make the bottles easily PE-recyclable.

Easy bottle handling

The screw caps have built-in cutting edges to open the bottles very easily. Open the cap, place it upside down on the bottle, press and turn the cap clockwise to break the ECO seal.

Easy automatic mixing

You can easily mix the parts using the automatic mixing units. Remove the cap, invert the bottle with the ECO seal over the seal opener and press. The seal is automatically broken and the contents can flow into the tank on the mixing unit. Whatever the mixing procedure, automatic or manual, mixing will never result in sludge formation.

Technical

SPECIFICATIONS

General rules

- Temperature adjustment of the developer by 1 °C equals speed adjustment by 2 to 3 log E points.

Suited for

- All kinds of processors, except Tabletop processors
- Not suited for manual processing.

Processing cycles

- All processing cycles, except HT

Replenishment rate

- Ranging from 300 to 600 ml/m²

Packaging

- One shipping carton E.O.S. Dev. contains developer for 2 x 20 l (2 x 5,28 gal)
- One shipping carton E.O.S. Fix contains fixer for 2 x 25 l (2 x 6,6 gal)

- One shipping carton E.O.S. Fixer contains fixer for 2 x 20 l (2 x 5,28 gal)

Conditions

- Make sure the processing room is well ventilated: air ventilation of 10 times the room volume per hour will prevent build-up of chemicals in the air.

Safety

- Consult the Safety Data Sheet (SDS) for information relating to safety, health and environmental aspects. Please follow the safety instructions mentioned on the bottles carefully.

Agfa and the Agfa rhombus are trademarks of Agfa-Gevaert N.V., Belgium, or its affiliates. All other trademarks are held by their respective owners and are used in an editorial fashion with no intention of infringement. The data in this publication are for illustration purposes only and do not necessarily represent standards or specifications, which must be met by Agfa HealthCare. All information contained herein is intended for guidance purposes only, and characteristics of the products and services described in this publication can be changed at any time without notice. Products and services may not be available for your local area. Please contact your local sales representative at agfa.com for availability information. Agfa HealthCare diligently strives to provide as accurate information as possible, but shall not be responsible for any typographical error.

© Copyright 2009 Agfa HealthCare NV

All rights reserved

Printed in Belgium

Published by Agfa HealthCare NV

B-2640 Mortsel-Belgium

NFWNL GB 00200908