



PHOTOCHEMISCHE FABRIK GMBH

CAWO X-Ray Protection

CAWO protective devices provide appropriate protection for staff and patients of radiographic departments. The aprons and adapted accessories are comfortable to wear and provide indispensable protection against X-rays.

Made for day-to-day use CAWO protective aprons consist of lead rubber layers covered by robust and non-irritant Nylon. Easy to handle and to clean. Regular quality-control and careful treatment guarantee security and a long life of the CAWO X-ray protection aprons.

CAWO-Flex ST: standard lead rubber
CAWO-Flex Light: newly developed interior material, up to 20 % less in weight

The standard colour for CAWO X-ray protection clothes is royal blue. We may offer a wide range of other colours on request.




 grey	 royal blue
 navy	 scarlet
 yellow	 petrol
 pink	 bordeaux
 purpur	 jade
 neon green	 orange

The CAWO protective devices comply with the requirements of the Council Directive 89/866/EEC concerning personal protective devices and are in conformity with the Standard IEC 61331-3. All aprons carry the CE mark and are delivered with a coat hanger.



PHOTOCHEMISCHE FABRIK GMBH





CAWO Staff Protection

<p>CAWO Coat Apron</p> <p>Long front part, short back part.</p> <p>Lead equivalents: 0.25 mm Pb front and back 0.35 mm Pb front, 0.25 mm Pb back 0.50 mm Pb front, 0.25 mm Pb back</p>		<p>Available sizes:</p> <table> <tbody> <tr> <td rowspan="3">Small</td> <td>60 x 90 cm</td> </tr> <tr> <td>60 x 100 cm</td> </tr> <tr> <td>60 x 110 cm</td> </tr> <tr> <td rowspan="3">Medium</td> <td>60 x 100 cm</td> </tr> <tr> <td>60 x 110 cm</td> </tr> <tr> <td>60 x 120 cm</td> </tr> <tr> <td rowspan="3">Large</td> <td>60 x 130 cm</td> </tr> <tr> <td>75 x 100 cm</td> </tr> <tr> <td>75 x 110 cm</td> </tr> <tr> <td></td> <td>75 x 120 cm</td> </tr> <tr> <td></td> <td>75 x 130 cm</td> </tr> </tbody> </table>	Small	60 x 90 cm	60 x 100 cm	60 x 110 cm	Medium	60 x 100 cm	60 x 110 cm	60 x 120 cm	Large	60 x 130 cm	75 x 100 cm	75 x 110 cm		75 x 120 cm		75 x 130 cm																				
Small	60 x 90 cm																																					
	60 x 100 cm																																					
	60 x 110 cm																																					
Medium	60 x 100 cm																																					
	60 x 110 cm																																					
	60 x 120 cm																																					
Large	60 x 130 cm																																					
	75 x 100 cm																																					
	75 x 110 cm																																					
	75 x 120 cm																																					
	75 x 130 cm																																					
<p>CAWO Double-sided Apron</p> <p>Large front and large back part.</p> <p>Lead equivalents: 0.25 mm Pb front and back 0.35 mm Pb front, 0.25 mm Pb back 0.50 mm Pb front, 0.25 mm Pb back</p>		<p>Available sizes:</p> <table> <tbody> <tr> <td>60 x 90 cm</td> </tr> <tr> <td>60 x 100 cm</td> </tr> <tr> <td>60 x 110 cm</td> </tr> <tr> <td>60 x 120 cm</td> </tr> </tbody> </table>	60 x 90 cm	60 x 100 cm	60 x 110 cm	60 x 120 cm																																
60 x 90 cm																																						
60 x 100 cm																																						
60 x 110 cm																																						
60 x 120 cm																																						
<p>CAWO Skirt and Jacket</p> <p>Combining two protective part for the upper and the lower part of the body a spread of weight to shoulders and hip is achieved. More comfort for the user.</p> <p>Lead equivalents: 0.25 mm Pb – 0.50 mm Pb in the front through overlapping</p>		<p>Skirt</p> <p>Available sizes:</p> <table> <thead> <tr> <th></th> <th>Length</th> <th>Width</th> </tr> </thead> <tbody> <tr> <td rowspan="3">Small</td> <td>55 cm</td> <td>65 – 75 cm</td> </tr> <tr> <td>60 cm</td> <td>65 – 75 cm</td> </tr> <tr> <td>65 cm</td> <td>65 – 75 cm</td> </tr> <tr> <td rowspan="3">Medium</td> <td>60 cm</td> <td>75 – 85 cm</td> </tr> <tr> <td>65 cm</td> <td>75 – 85 cm</td> </tr> <tr> <td>70 cm</td> <td>75 – 85 cm</td> </tr> <tr> <td rowspan="3">Large</td> <td>60 cm</td> <td>85 – 95 cm</td> </tr> <tr> <td>65 cm</td> <td>85 – 95 cm</td> </tr> <tr> <td>70 cm</td> <td>85 – 95 cm</td> </tr> </tbody> </table> <p>Jacket</p> <p>Available sizes:</p> <table> <thead> <tr> <th></th> <th>Length</th> <th>Chest</th> </tr> </thead> <tbody> <tr> <td>Small</td> <td>60 cm</td> <td>107 cm</td> </tr> <tr> <td>Medium</td> <td>60 cm</td> <td>117 cm</td> </tr> <tr> <td>Large</td> <td>60 cm</td> <td>127 cm</td> </tr> </tbody> </table>		Length	Width	Small	55 cm	65 – 75 cm	60 cm	65 – 75 cm	65 cm	65 – 75 cm	Medium	60 cm	75 – 85 cm	65 cm	75 – 85 cm	70 cm	75 – 85 cm	Large	60 cm	85 – 95 cm	65 cm	85 – 95 cm	70 cm	85 – 95 cm		Length	Chest	Small	60 cm	107 cm	Medium	60 cm	117 cm	Large	60 cm	127 cm
	Length	Width																																				
Small	55 cm	65 – 75 cm																																				
	60 cm	65 – 75 cm																																				
	65 cm	65 – 75 cm																																				
Medium	60 cm	75 – 85 cm																																				
	65 cm	75 – 85 cm																																				
	70 cm	75 – 85 cm																																				
Large	60 cm	85 – 95 cm																																				
	65 cm	85 – 95 cm																																				
	70 cm	85 – 95 cm																																				
	Length	Chest																																				
Small	60 cm	107 cm																																				
Medium	60 cm	117 cm																																				
Large	60 cm	127 cm																																				



PHOTOCHEMISCHE FABRIK GMBH

CAWO Patient Protection

<p>CAWO Protective Gonad Apron</p> <p>The CAWO-Flex protective gonad apron is for the protection of patients during chest radiography or screening.</p> <p>Lead equivalence: 0.50 mm Pb</p>		<p>Available sizes: 20 x 25 cm 30 x 30 cm 37 x 45 cm 40 x 50 cm</p>
<p>CAWO Scrotum Shield</p> <p>Lead equivalence: 1.0 mm Pb</p>		<p>Available sizes: small medium large</p>
<p>CAWO Ovary Shield</p> <p>Lead equivalence: 1.0 mm Pb</p>		<p>Set of 3 sizes (small, medium, large)</p>
<p>CAWO Protective Thyroid Collar</p> <p>Lead equivalence: 0.50 mm Pb</p>		<p>Available sizes: Standard or with extended front protection</p>
<p>CAWO Dental Apron</p> <p>Lead equivalence: 0.25 mm Pb front and back 0.35 mm Pb front and back 0.50 mm Pb front and back</p>		<p>Available sizes: 60 x 80 cm</p>



PHOTOCHEMISCHE FABRIK GMBH

CAWO Patient Protection

CAWO Image Intensifying Apron

0.25 / 0.35 mm Pb



Available sizes:

60 x 100 cm

60 x 110 cm



60 x 120 cm

60 x 130 cm



PHOTOCHEMISCHE FABRIK GMBH

CAWO Protective Accessories

<p>Protective Gloves with inside gloves</p> <p>Lead equivalence: 0.25 mm Pb, grey 0.50 mm Pb, green</p>		<p>One size for all</p>
<p>Lead Glass Spectacles</p> <p>Lead equivalence: 0.75 mm Pb</p>		<p>One size for all</p>
<p>Protective Cover Sheet (lead rubber cassette covering)</p> <p>Lead equivalence: 0.80 mm Pb</p>		<p>Sizes as required</p>
<p>Lead Glass for laboratories</p> <p>Lead equivalence: 1.6 mm Pb 2.3 mm Pb 2.8 mm Pb 3.3 mm Pb 3.6 mm Pb</p>		<p>Available sizes: Sizes as required but max. 100 x 200 cm</p>

CAWO radiographic protection accessories comply with the requirements of the Council Directive 89/686/EEC concerning personal protective devices and are in conformity with the Standard IEC51331-3. They carry the CE-mark.